## »VPS Putty«

de

### Produktbeschreibung

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en-US

### Product description

*VPS Putty/VPS Putty (fast) VPS Putty (soft)* is a VPS precision impression material with a high consistency and is based on addition-curing silicones.

### Indications

Dental cases that require the production of dental models.

### Contraindications

Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended use

Impressions using the correction impression or double-mixing technique combined with *VPS Light* for the fabrication of crowns, bridges, inlays, and onlays.

### Limitations of use

Use the material for two-step impressions only.

### Intended users

Dentist, dental assistant staff

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Recommended use

#### Preparing the impression

Have *VPS Light* correction impression material ready for use and, shortly before application, check that the openings of the safety cartridges are not obstructed. The used mixing tip acts as a seal until the next application is carried out. Do not reuse the original cap of the safety cartridge under any circumstances! This can result in unwanted curing due to contamination of the base and catalyst pastes. If retraction fluids have been used, carefully rinse and dry the region prior to taking the impression.

#### Mixing the materials

Mix at a ratio of 1:1 and knead with your fingertips for at least 30 seconds until the color is homogeneous. Do not use latex gloves. Close containers with the correct lids immediately after use. Do not mix up the container lids.

#### Corrective impression

Insert the impression tray filled with *VPS Putty/VPS Putty (fast)/VPS Putty (soft)* into the mouth. This position should be maintained for at least as long as the specified setting time in the mouth for the material. After removal, thoroughly clean and dry the initial impression. Prior to carrying out the final impression, cut relief channels into the initial impression and eliminate any imperfections.

Fill *VPS Light* correction impression material into the initial impression, and if necessary, inject the material around the prepared teeth. Place the filled initial impression back into the mouth, observing the working and setting time of the correction impression material. To ensure that the two impression materials bond securely together, do not wait longer than 30 min between the completion of the initial impression and the preparation of the final impression.

#### Double-mixing technique

While the assistant loads the impression tray, the dentist can inject *VPS Light* around the prepared teeth. Try to complete the filling of the impression tray and the injection around the prepared teeth at the same time, since the higher temperature in the patient's mouth causes the correction impression material to cure more quickly than the tray material. Place the impression tray into the mouth. This position should be maintained for at least as long as the specified setting time in the mouth for the material.

#### Disinfection and casting

The impression can be disinfected using commercially available dipping solutions. The use of a 2% glutaraldehyde solution is recommended. Please ensure that the disinfectant used is indicated for use with addition-cured silicone impression materials.

The exposure period for the impression in the disinfectant is approx. 10 - 15 min. Afterwards, rinse the impression under running water for 15 s. Store the impression at normal room temperature (max. 25 °C/77°F) and wait at least 30 min after removing the impression from the mouth.\*

The impression material is compatible with all commercially available plasters. The use of a hard plaster (ISO 6873:2013, Type 3, e.g. Moldano, Bayer) for the cast and an ultra-hard dental plaster for the molds (ISO 6873:2013, Type 4, e.g. Fujirock II, GC) is recommended. The impressions can be copper or silver-plated.

\* Additional information: Statement by the DGZMK 3/93 V 1.0; yiry.dgzmk.de ADA recommends following the advice of the Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA1991;122(3):110; wym.ada.org

### Timing

VPS Putty	
Mixing time	30 s
Working time (from the start of mixing)*	≤ 105 s
Recommended setting time in the mouth	≥ 225 s

VPS Putty (fast)	
Mixing time	30 s
Working time (from the start of mixing)*	≤ 75 s
Recommended setting time in the mouth	≥ 180 s

VPS Putty (soft)	
Mixing time	30 s
Working time (from the start of mixing)*	≤ 90 s
Recommended setting time in the mouth	≥ 225 s

\* The specified working times apply at a room temperature of 23 °C (72 °F) and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.  
The working time describes the period after mixing has begun until the time by which the impression tray must be positioned in the patient's mouth.

## »VPS Putty«

### Notes for use

- The recommended times can be found in the section headed "Timing".
- Perforated or non-perforated metal or resin trays can be used with the impression material.
- When selecting the impression tray, ensure that there is enough room between the wall of the impression tray and the dental arch.
- In order to improve the adhesion of the impression material to non-perforated
- impression trays, apply a thin, even film of *Tray adhesive for VPS*. Allow *Tray adhesive for VPS* to dry for 3 min before loading the impression tray.

### Interactions

- Do not use with condensation-curing or polyether impression materials.
- Do not allow the impression to come into contact with organic solvents or liquids that contain solvents as this can cause the material to swell and lead to inaccuracies.
- As retraction solutions may impair the setting time of the material, any retraction solution residues must be removed prior to taking the impression by thoroughly spraying and drying the region.
- Do not use with latex gloves because they inhibit the setting process. The use of PVC or PE-based gloves is recommended.

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings/precautions

- Keep out of the reach of children!
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Product performance characteristics

The product complies with the requirements of ISO 4823:2021: putty consistency.

### Storage/disposal

- Store in a dry place at 15 to 25°C (59 to 77°F)!
- Protect against exposure to direct sunlight!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

### Composition

REF 9043MDR: SiO<sub>2</sub>, vinyl silicones/siloxanes, methyl silicones/siloxanes, hydrogen silicone, paraffin/vaseline, additives, pigments, catalysts.

REF 9044MDR: SiO<sub>2</sub>, vinyl silicones/siloxanes, methyl silicones/siloxanes, additives, pigments, catalysts, hydrogen silicone.

REF 9045MDR: SiO<sub>2</sub>, vinyl silicones/siloxanes, methyl silicones/siloxanes, additives, pigments, catalysts, hydrogen silicone.

### Packaging

<b>VPS Putty</b>	
REF 9043MDR	2 Jars @ 262 ml
<b>VPS Putty (fast)</b>	
REF 9044MDR	2 Jars @ 262 ml
<b>VPS Putty (soft)</b>	
REF 9045MDR	2 Jars @ 262 ml

**»VPS Putty«**