»VPS Light/VPS Light (fast)«

de

Produktbeschreibung

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Das Produkt entspricht ISO 4823:2021, low consistency-light bodied

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en-US

Product description

VPS Light/VPS Light (fast) is a low-consistency precision impression material based on VPS addition-cured silicone and is automatically mixed.

Indications

Dental cases that require the production of dental models.

Contraindications

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended use

Impressions using correction impression or double-mixing technique combined with *VPS Putty* for fabricating crowns, bridges, inlays and onlays.

Limitations of use

Use the material for two-phase impressions only.

Intended users

Dentist, dental assistant staff

Patient target group

Persons receiving treatment as part of a dental procedure.

Recommended use

Preparing the impression

Have *VPS Light/VPS Light (fast)* correction impression material ready for use, and shortly before application, check that the openings of the safety cartridge are not obstructed. The used mixing tip acts as a seal until the next application is carried out. Do not reuse the original closure cap of the safety cartridge under any circumstances! This can result in unwanted curing due to contamination of the base and catalyst pastes. If retraction fluids have been used, carefully rinse and dry the region prior to taking the impression.

Corrective impression

Load the impression tray with *VPS Putty* and place the tray into the patient's mouth. This position should be maintained for at least as long as the specified setting time in the mouth for the material. After removal, thoroughly clean and dry the initial impression. Prior to performing the correction impression, cut drainage channels into the initial impression and eliminate any defects.

Place *VPS Light/VPS Light (fast)* correction impression material into the initial impression and if necessary, inject the material around the prepared teeth. Place the filled correction impression back into the mouth, observing the working and setting time of the correction impression material. To ensure that the two impression materials bond securely together, do not wait longer than 30 min between the completion of the initial impression and preparation of the final impression.

Double-mixing technique

While the assistant loads the tray, the dentist can inject around the prepared teeth with *VPS Light/VPS Light (fast)*. Try to complete the filling of the impression tray and the injection around the prepared teeth at the same time, since

the higher temperature of the patient's mouth causes the correction impression material to cure more quickly than the tray material. Place the impression tray into the mouth. This position should be maintained for at least as long as the specified setting time in the mouth for the material.

Disinfection and casting

The impression can be disinfected using commercially available dipping solutions. The use of a 2 % glutaraldehyde solution is recommended. Please ensure that the disinfectant used is indicated for use with addition-cured silicone impression materials.

The exposure period for the impression is approx. 10-15 min. Afterwards, rinse the impression under running water for 15 s. Store the impression at normal room temperature (max. 25 °C/77 °F). Wait at least 30 min after removing the impression from the patient's mouth before casting.*

The impression material is compatible with all commercially available plasters. The use of a hard plaster (ISO 6873:2013, Type 3, e.g. Moldano, Bayer) for the cast and an ultra-hard dental plaster for the molds (ISO 6873:2013, Type 4, e.g. Fujirock II, GC) is recommended. The impressions can be copper or silver-plated.

* For more information, see the statement issued by the German Society for Dental and Oral Medicine (DGZMK). Comments of the DGZMK 3/93 V 1.0; www.dgzmk.de ADA recommends following the advice of the Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA 1991;122(3):110; www.ada.org

Storing the impression

- Store the impression at 15–25 °C (59–77 °F) in a dry place!
- Protect against exposure to direct sunlight!
- The impression is stable and can be cast even after months without loss of accuracy.

Timing

VPS Light	
Working time (from the start of mixing)*	≤ 120 s
Recommended setting time in the mouth	≥ 210 s
VPS Light (fast)	
Working time (from the start of mixing)*	≤ 90 s
Recommended setting time in the mouth	≥ 135 s

* The specified working times are valid for a room temperature of 23 °C (72°F) and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

The working time describes the period after mixing has begun until the time by which the impression tray must be positioned in the patient's mouth.

Notes for use

- The recommended times can be found in the section headed »Timing«.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Interactions

- Do not use with condensation-curing or polyether impression materials.

»VPS Light/VPS Light (fast)«

- Do not allow the impression to come into contact with organic solvents or liquids that contain solvents as this can cause the material to swell and lead to inaccuracies.
- As retraction solutions may impair the setting time of the material, any retraction solution residues must be removed prior to taking the impression by thoroughly spraying and drying the region.
- Do not use with latex gloves because they inhibit the setting process. The use of PVC or PE-based gloves is recommended.

Cartridge handling

The double syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new cartridge, remove the transport cap first. The mixing tip and, if necessary, an intraoral tip can then be fitted. After use and until the next treatment, leave the mixing tip on the cartridge as a seal. The used mixing tip is replaced with a new tip when the device is used for the next treatment.

Note: When using a double syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g. a pea-sized amount).

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

The product complies with ISO 4823:2021, low consistency, light-bodied

Storage/disposal

- Store in a dry place at 15 to 25°C (59 to 77°F).
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

REF 9054MDR: Vinyl silicones/siloxanes, methyl silicones/siloxanes, SiO₂, hydrogen silicone, additives, pigments, catalysts.

REF 9055MDR: Vinyl silicones/siloxanes, methyl silicones/siloxanes, SiO₂, hydrogen silicone, additives, pigments, catalysts.

Packaging

VPS Light	
REF 9054MDR	2 Cartridges @ 50 ml 12 Mixing Tips 12 Intraoral-Tips
VPS Light (fast)	
REF 9055MDR	2 Cartridges @ 50 ml 12 Mixing Tips 12 Intraoral-Tips