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## Instructions for use » VPS Putty (fast)«

VPS Putty (fast) is a VPS precision impression material with high consistency, based on addition-cured silicone.

### Indication

Correction impression or double mixing technique in combination with VPS Light for fabricating crowns, bridges, inlays and onlays.

### Impression trays and adhesive

Perforated or non-perforated, metal or resin trays can be used with VPS Putty (fast). When selecting an impression tray, ensure that there is adequate space available between the tray wall and dental arch. Tray adhesive for VPS should be applied as a thin, uniform film to improve the adhesion of VPS Putty (fast) to non-perforated trays. Allow the Tray adhesive for VPS to dry for 3 minutes before loading the impression tray with material.

### Cartridge handling

Place the cartridge in the mixing device. When doing so observe the instructions for use issued by the respective device manufacturer.

### Attaching / changing the mixing tip

1. Open the clip by sliding it along its guide until it reaches the stop point and then drawing the original seal or the old mixing tip off. (Figure 1-2)


Figure 1


Figure 2

2. Now set a new mixing tip in place and press both openings on the rear side of the mixing tip firmly over the respective outlets on the cartridge. (Figure 3)

Figure 3

3. Lock the mixing tip in place with the aid of the clip, by sliding it along its guide until it encloses the mixing tip and reaches the stop position. (Figure 4)

Figure 4

Note: Always attach a new mixing tip prior to every mixing process! The mixing devices are not able to squeeze the impression material through cured material remaining in the mixing tip.

Note: Ensure that the mixing tip is seated firmly on the cartridge! Otherwise the impression material will be squeezed past the mixing tip, which may produce messy results.

### Recommended use

#### Preparing the impression

Prepare the VPS Light correction material and check that the outlets of the safety cartridges are free from obstruction shortly before application. Used mixing tips serve as a seal until the next application. Do not use the original seal on the safety cartridge again! Multiple use can cause unintentional curing due to contamination of the base and catalyst paste. If retraction liquids have been used, thoroughly rinse and dry the treatment site prior to taking the impression.

#### Correction impression

Insert the impression tray loaded with VPS Putty (fast) into the patient's mouth. Maintain this position until the recommended time in mouth of the material. Thoroughly clean and dry the impression after removal. Prepare escape channels in the preliminary impression and remove high spots before taking the correction impression.

Apply VPS Light correction material to the pre-impression and syringe correction material around the prepared teeth if necessary. Re-insert the loaded preliminary impression into the patient's mouth observing the working and setting times of the correction material. Do not wait longer than 30 minutes between the preliminary impression and correction impression to ensure that the two impression materials bond together.

#### Double mixing technique

While the assistant loads the tray, the operator can syringe around the prepared teeth with VPS Light. Loading the tray and syringing around the prepared teeth should be completed simultaneously if possible, as the correction material sets quicker than the tray material due to the higher intraoral temperature. Insert the impression tray into the patient's mouth. Maintain this position for the recommended time in mouth of the material.

#### Disinfecting and casting the impression

The impression can be disinfected with commercially available immersion solutions. We recommend using a 2% glutaraldehyde solution. Ensure that the disinfectant is indicated for use with addition-cured impression silicones.

Immerse the impression for approx. 10 to 15 minutes. Then rinse under running water for 15 seconds. Keep the impression at normal room temperature (max. 25 °C / 77 °F) and cast the impression no sooner than 30 minutes after removing it from the patient's mouth.\*

The impression material is compatible with all commercial casts. We recommend using a special purpose hard plaster (ISO 6873, Type 3; e.g. Moldano, Bayer) for the cast or an ultra hard dental plaster for the forms (ISO 6873, Type 4; e.g. Fujirock II, GC). Impressions can be copper-plated or silver-plated.

\* For further information see the comments of the DGZMK. Comments of the DGZMK 3/93 V 1.0; [www.dgzmk.de](http://www.dgzmk.de)

ADA recommends following Council on Dental Materials, Instruments and Equipment. Disinfection of impressions. JADA 1991;122(3):110; [www.ada.org](http://www.ada.org)

#### Please note

- Do not use with condensation-cured or polyether impression materials.
- Avoid contact of the impression with organic solvents or liquids containing solvent, as this can cause the material to swell and lead to inaccuracies.
- Retraction solutions should be thoroughly removed through spraying and drying before taking the impression, to avoid them impairing the setting properties. Do not use with latex gloves.
- In rare cases it cannot be ruled out that predisposed patients may develop hypersensitivity to individual components of the material.

### **Technical data**

#### **VPS Putty (fast)**

Working time (incl. mixing time)  $\leq$  1:45 min.  
Recommended time in mouth  $\geq$  3:15 min.  
Strain in compression approx. 1.3 – 2.3 %  
Compression set approx. 0.7 %  
Shore-A hardness (after 24 h) approx. 75

Times apply at a normal room temperature of 23 °C / 73 °F and a normal relative humidity of 50%. Higher temperatures shorten, lower temperatures increase those times.

#### **Composition**

Vinyl polysiloxanes (VPS), silicon dioxide, platinum catalyst, paraffin, tensides and small quantities of stabilizers and additives.

#### **Classification**

DIN EN ISO 4823:2000, type 0, putty consistency

#### **Storage**

Store at room temperature (15 - 25 °C / 59 - 77 °F). Do not use after the expiration date.

#### **Available product sizes**

##### **VPS Putty (fast)**

REF 9068 1 cartridge @ 380 ml

**Keep out of reach of children! For dental use only!**