Gebrauchsinformation »Nano hybrid composite«

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Handelsformen

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REF 2061 20 Composite Tips (A2) à 0,3 g

Instructions for use for "Nano hybrid composite"

Nano hybrid composite is a light-cured, radiopaque filling material for restorations of all cavity classes.

Indications

Restorations of all cavity classes

Contraindications

- Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.
- Do not apply directly to exposed pulp or use in cases of known allergy to any of the components.

Interactions with other materials

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be excluded in individual cases.

Basic safety information

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

Instructions for use

- Light units should have an output of 450 nm and should be checked regularly. The light intensity should be at least 400 mW/cm². Place the light as close as possible to the filling material.
- Do not remove the oxygen inhibited layer, because this is required for the bond with the next layer.
- Follow the manufacturer's instructions for any other materials that are used with *Nano hybrid composite*.

Recommended use

1. Clean the tooth. Remove all cleaning residue with water spray. Determine the VITA shade. The use of a rubber dam is recommended.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
3. If necessary, etch using your preferred etching technique such as Total-Etch. Observe the manufacturer's instructions for the material used.
4. Apply bonding agent in accordance with the manufacturer's instructions.
5. If necessary, apply an underfilling of low-viscosity composite or compomer. The use of a matrix is recommended.
6. Apply the desired shade of *Nano hybrid composite* to the cavity with a maximum layer thickness of 2 mm. Light-cure each surface for 20 seconds
7. Removal of any excess material, contouring, and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size. *Nano hybrid composite* can be given a high shine using polishing paste together with polishing cups or discs.

Composition

Barium glass in a Bis-GMA-based matrix of dental resins. Pigments, additives, and catalysts. Filler content: 82% by weight = 65% by volume. The range of variation of the inorganic filler particles is between 0.02-3.0 µm.

Classification

DIN EN ISO 4049:2009, Type 1, Class 2

Storage

Store at room temperature (15-25 °C/59-77 °F). Close syringes/composite-Tip immediately after use so that the material does not cure in the ambient light. Do not use after the expiration date.

Available product sizes

REF 2054 1 syringe (A2) containing 3 g

REF 2061 20 composite tips (A2) containing 0.3 g