»Universal composite luting cement«

de

Produktbeschreibung

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en-us

Product description

Universal composite luting cement is a universally applicable, dual-curing, bonding composite luting cement for the permanent insertion of indirect restorations and root posts. The radiopaque material is used together with a suitable bonding agent (e.g. *Universal adhesive dual*). *Universal composite luting cement* is automatically mixed in the double-barreled syringe and applied directly to the restoration or into the cavity. Manual mixing of the two components is not necessary.

Indications

- Extensive tooth defects: Tooth decay, severe wear, traumatic injuries of the dental enamel
- Internal tooth discoloration, dental misalignments and deformities
- Tooth loss

Contraindications

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended use

- Permanent luting of crowns, inlays, onlays and bridges from metal and metal alloys (precious and base metals), oxide ceramics (e.g. zircon) and all silicate ceramics (e.g. lithium disilicate, feldspar) as well as composites.
- Permanent cementation of ceramic or composite veneers
- Permanent luting of oxide ceramic, silicate ceramic, composite or metal restorations on implant abutments
- Cementation of root posts

Limitations of use

- Do not use the material for direct capping of the pulp material.
- Do not use the material if a sufficiently dry working area or the recommended application technique is not possible

Intended users

Dentist

Patient target group

Persons receiving treatment as part of a dental procedure.

Recommended use

Preparing the tooth substance

1. Prepare the tooth substance in accordance with the general rules. Prepare the root canal endodontically and prepare for the root post insertion. Leave at least 4 mm gutta-percha apically.
2. Remove any residual temporary luting material if necessary.
3. Afterwards, rinse the prepared tooth and dry it with air. Avoid both excessive residual moisture and over-drying.
4. Try in the restoration

Pre-treating the restoration materials

Oxide ceramic/metal

Please observe the following: Do not etch oxide ceramic.

1. Pre-treat the restoration in accordance with the manufacturer's instructions. If not otherwise described, sandblast the surfaces to be bonded with aluminum oxide ($\leq 50 \mu\text{m}$).
2. Clean with alcohol and dry with air that is free of oil and water.
3. Pre-treat using a suitable primer in accordance with the manufacturer's instructions.

Etchable silicate ceramic

1. Pre-treat the restoration in accordance with the manufacturer's instructions. If not otherwise described, rinse the surfaces to be bonded thoroughly with 5% liquid acid as per manufacturer's instructions.
2. Dry with air that is free of oil and water.
3. Prepare the surfaces to be bonded with a silanizing agent as per manufacturer's instructions.

Composite restorations

1. Pre-treat the restoration in accordance with the manufacturer's instructions. If not described elsewhere, sandblast the surfaces to be bonded with aluminum oxide until they have a matt appearance.
2. Clean with alcohol and dry with air that is free of oil and water.
3. Apply the primer.

Glass-fiber reinforced root post

Pre-treat in accordance with the manufacturer's instructions.

Apply cement to the pre-treated tooth substance or restoration

Please observe the following: Use of a rubber dam is recommended for isolating the working area.

1. Apply a suitable adhesive, e.g., *Universal adhesive dual*, in accordance with the manufacturer's instructions.
2. Add *Universal composite luting cement* to the pre-treated restoration

Please observe the following: Given the short curing time it is advisable to cement a maximum of two veneers at a time.

3. Insert restoration within 01:00 min after start of mixing and begin luting onto the prepared tooth with light pressure and stabilize until the cement has cured.
4. Remove excess within 00:30 to 01:00 min of inserting the restoration.
5. Remove interdental excess carefully using dental floss. Alternatively: Briefly light-cure any excess cement (1-2 s) and remove.
6. Once mixing begins, allow the material to chemically cure completely for 07:00 min.
7. In the case of light-permeable restorations, additionally cure the material with a suitable light unit for at least 20 seconds.

»Universal composite luting cement«

Apply cement to the pre-treated root canal

Please observe the following: Use of a rubber dam is recommended for isolating the working area.

Please observe the following: The working time of *Universal composite luting cement* can be reduced by contact with a dual-curing adhesive.

1. Apply a suitable adhesive, ideally a dual-curing one, as per manufacturer's instructions.
2. Apply *Universal composite luting cement* directly to the pre-treated root canal using a suitable tip. Introduce the tip as deeply into the root canal as possible.

Please observe the following: Leave the tip immersed in the material and do not pull up until the material starts to rise.

3. Insert the root post within 00:30 min after mixing has begun and cement using slight pressure.
4. Remove excess within 00:30 min of inserting the root post.
5. Once mixing begins, allow the material to chemically cure completely for 07:00 min.

Note: The material can be additionally cured with a suitable light unit for at least 20 seconds.

Notes for use

- Light units should have a light intensity of 600 mW/cm², emit at 450 nm, and their performance should be checked regularly. The light-curing unit should be positioned as closely to the material as possible. Observe the manufacturer's instructions.
- It is advisable to use a rubber dam.
- It is important that the surface to be treated is not dirty before each following working step.
- During polymerization of the cement, an oxygen inhibition layer develops in the areas that have been exposed to atmospheric oxygen. This layer does not set. For this reason, only large areas of excess material should be removed, allowing smaller areas of excess material to remain. The oxygen inhibition layer is removed together with the smaller areas of excess material when finishing the edges.
- Observe the instructions for use of the materials that are used with the *Universal composite luting cement*.

Interactions

- Water and air containing oil can prevent the polymerization of the materials at the contact surfaces.
- Not to be used in conjunction with materials that contain eugenol / clove oil as these inhibit the curing of the materials.
- Hydrogen peroxide can impair the curing of the material.
- Phosphoric acid can hinder an optimal bond with zircon oxide aluminum oxide, and metal restorations.
- The working time of the luting cement can be reduced by contact with an adhesive.

Using the double-barreled syringe

The double-barreled syringe enables automatic dispensing and mixing of the material, which allows a time-saving direct application. When using a new double-barreled syringe, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the double-barreled syringe as a seal. The used mixing tip is replaced with a new tip when the device is used for the next treatment.

Note: When using a double-barreled syringe for the first time, discard a small amount of the initially ejected material (e.g. a pea-sized amount).

Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Performance characteristics of the device

The product complies with ISO 4049: 2019, Type 2, Class 3 (incl. shade stability).

Storage / disposal

- Store in a dry place in the fridge at a temperature of 2 – 25 °C/36 – 77 °F
- Protect against exposure to direct sunlight!
- After application, leave the used mixing tip on the double-barreled syringe as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

Dental glass, EBPDMA, UDMA, TEDMA, HEMA, SiO₂, bis-GMA, mmHE, additives. Inorganic filling materials: approx. 51 vol% (0.02 -3 µm).

Packaging

| | |
|-------------|---|
| REF 2108MDR | 1 Double-barreled syringe (A2.5) @ 5 ml paste Accessories |
|-------------|---|