Gebrauchsinformation »Universal composite luting cement«

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Instructions for use »*Universal composite luting cement*«

Universal composite luting cement is a universally applicable, dual-curing, bonding composite luting cement for the permanent insertion of indirect restorations and root posts. The radiopaque material is used together with a suitable bonding agent (e.g. *universal adhesive dual*). *Universal composite luting cement* is automatically mixed in the double-barreled syringe and applied directly to the restoration or into the cavity. Manual mixing of the two components is not necessary.

Indications

- Permanent luting of crowns, inlays, onlays and bridges from metal and metal alloys (precious and base metals), oxide ceramics (e.g. zircon) and all silicate ceramics (e.g. lithium disilicate, feldspar) as well as composites
- Permanent luting of ceramic or composite veneers
- Permanent luting of oxide ceramic, silicate ceramic, composite or metal restorations on implant abutments
- Cementation of root posts

Contraindications

- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use the material if a sufficiently dry working area or the recommended application technique is not possible.
- Do not use the material for direct capping of the pulp material.

Basic safety information

- For dental use only.
- Keep out of the reach of children.
- Avoid skin contact with the paste. In case of accidental contact, rinse with water and soap. In case of contact with the eyes, rinse immediately with lots of water and consult a physician.

Side effects

There are no known systemic side effects to date. Rare cases of hypersensitivity to components of the material cannot be excluded. Do not use the material for direct capping of the pulp material.

Interactions

- Water and air containing oil can prevent the polymerization of the materials at the contact surfaces.
- Not be used in conjunction with materials that contain eugenol / clove oil as these inhibit the hardening of the materials.
- Hydrogen peroxide can impair the hardening of the material.
- Phosphoric acid can compromise the optimal bond with zircon oxide aluminum oxide and metal restorations.
- The working time of the luting cement can be reduced by contact with an adhesive.

Notes for use

- Light units should have a luminous intensity of 400 mW/cm², emit at 450 nm, and their performance should be checked regularly. The light should be positioned as closely to the material as possible. Observe the manufacturer's instructions.
- The use of a rubber dam is recommended.
- It is important that the surface to be treated is not dirty before each following working step.
- During polymerization of the cement, an oxygen inhibition layer develops in the areas that have been exposed to atmospheric oxygen. This layer does not set. For this reason, only large areas of excess material should be removed, allowing smaller areas of excess material to remain. The oxygen inhibition layer is removed together with the smaller areas of excess material when finishing the edges.
- Observe the instructions for use of the materials that are used with the *Universal composite luting cement*.

Recommended use

Preparing the tooth substance

1. Prepare the tooth substance in accordance with the general rules. Prepare the root canal endodontically and prepare for the root post insertion. Leave at least 4 mm gutta-percha apically.
2. Remove any residual temporary luting material if necessary.
3. Afterwards, rinse the prepared tooth and dry it with air. Avoid both excessive residual moisture and over-drying.
4. Try in the restoration.

Pre-treating the restoration materials

Oxide ceramic/metal

Please observe the following: Do not etch oxide ceramic.

1. Pre-treat the restoration in accordance with the manufacturer's instructions. If not otherwise described, sandblast the surfaces to be bonded with aluminum oxide (≤50 µm).
2. Clean with alcohol and dry with air that is free of oil and water.
3. Pre-treat using a suitable primer in accordance with the manufacturer's instructions.

Etchable silicate ceramic

1. Pre-treat the restoration in accordance with the manufacturer's instructions. If not otherwise described, rinse the surfaces to be bonded thoroughly with 5% liquid acid as per manufacturer's instructions.
2. Dry with air that is free of oil and water.
3. Prepare the surfaces to be bonded with a silanizing agent as per manufacturer's instructions.

Composite restorations

1. Pre-treat the restoration in accordance with the manufacturer's instructions. If not described elsewhere, sand-blast the surfaces to be bonded with aluminum oxide until they have a matt appearance.
2. Clean with alcohol and dry with air that is free of oil and water.
3. Apply the primer.

Glass-fiber reinforced root post

- Pre-treat in accordance with the manufacturer's instructions.

Apply cement to the pre-treated tooth substance or restoration

Please observe the following: Use of a rubber dam is recommended for isolating the working area.

1. Apply a suitable adhesive, e.g. *Universal adhesive dual*, in accordance with the manufacturer's instructions.
2. Add *Universal composite luting cement* to the pre-treated restoration.

Please observe the following: Given the short curing time it is advisable to cement a maximum of two veneers at a time.

3. Insert restoration within 1:00 min after start of mixing and begin luting onto the prepared tooth with light pressure and stabilize until the cement has cured.
4. Remove excess within 0:30 to 1:00 min of inserting the restoration.
5. Remove interdental excess carefully using dental floss. Alternatively: Briefly light-cure any excess cement (1-2 s) and remove.
6. Once mixing begins, allow the material to chemically cure completely for 7:00 min.
7. In the case of light-permeable restorations, additionally cure the material with a suitable light unit for at least 20 seconds.

Apply cement to the pre-treated root canal

Please observe the following: Use of a rubber dam is recommended for isolating the working area.

Please observe the following: The working time of the *Universal composite luting cement* can be reduced by contact with a dual-curing adhesive.

1. Apply a suitable adhesive, ideally a dual-curing one, as per manufacturer's instructions.
2. Apply *Universal composite luting cement* directly to the pre-treated root canal using a suitable tip. Introduce the tip as deeply into the root canal as possible.

Please observe the following: Leave the tip immersed in the material and do not pull up until the material starts to rise.

3. Insert the root post within 0:30 min after mixing has begun and start luting using slight pressure.
4. Remove excess within 0:30 min of inserting the root post.
5. Once mixing begins, allow the material to chemically cure completely for 7 min.

Note: The material can be additionally cured with a suitable light unit for at least 20 seconds.

Using the double-barreled syringe

The material is dispensed and mixed automatically by extruding through a mixing tip. When using a new double-barreled syringe, remove the cap. Then attach the mixing tip and a tip if necessary. When using a double-barreled syringe for the first time, squeeze out a pea-sized amount of material and discard it.

After use, the mixing tip acts as a cap and should be left on the double-barreled syringe until the next time it is used. When the double-barreled syringe is used again, the mixing tip will then need to be replaced with a new one. Do not reuse the original cap of the double-barreled syringe under any circumstances. This can result in unwanted curing due to contamination of the base and catalyst paste.

Storage

Store in a dry place in the fridge at a temperature of 2–25 °C/36–77 °F. Keep away from direct sunlight. Do not use the material after the expiration date. Close packaging immediately after use.

Composition

Barium glass in a Bis-GMA-based matrix of dental resins, pigments, additives and catalysts.
Filler content: 69% by weight = 51 by vol. % (0.02 – 3.0 µm)

Classification

ISO 4049:2009 (incl. shade stability)

Packaging

REF 2108 1 double-barreled syringe (A2.5) @ 5 ml
paste, accessories