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Instructions for use » *Varnish & Bond 2*«

Methacrylate-based *Varnish & Bond 2* is used for sealing the surface and repair of temporary crowns and bridges.

Indication

Sealing the surface and repair of temporary restorations as well as primer/ bonding agent for indirect restorations of PEEK, composite and PMMA.

Recommended application for surface sealing

1. Apply a thin film of *Varnish & Bond 2* to the surface to be sealed and leave to work for 20 s.
2. Cure the material:
 - Cure the material per bridge pontic and per crown according to the times given below (see technical data).
 - When using a light-curing unit, the curing times given below (see technical data) apply for the entire restoration.
3. In both cases, cure the material within 2:00 min. to prevent formation of an oxygen inhibition layer.

Recommended application when repairing temporary restorations

1. Roughen the surface to be repaired with a rotary cutter when using as a bonding agent.
2. Apply a thin film of *Varnish & Bond 2* to the prepared surface.
3. Light cure the material for 20 s per joint, using a halogen lamp or an LED lamp. Alternatively, leave the entire restoration to cure in the light-curing unit for a minimum 90 s.
4. Apply freshly mixed temporary material or light-curing composite to the surface to be repaired (adhere to the instructions for use of the respective material).

Bonding agent for restorations

Bonding agent for restorations

1. Prepare the restoration according to manufacturer's specifications.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces. Expose for 20 s with a halogen lamp or ≥ 10 s with an LED lamp (light intensity $> 1000 \text{ mW/cm}^2$).
1. Lute restoration with a suitable adhesive dual- or self-curing composite luting cement according to the manufacturer's instructions.

Please note

- Position the light as near as possible to the material.
- *Varnish & Bond 2* cannot be used as dentine or enamel adhesive.
- Adhere to the instructions for use for other materials used with *Varnish & Bond*.
- *Varnish & Bond 2* also bonds to teeth as well as polycarbonate and polyacetate foil.
- Carefully seal the bottle immediately after each use to prevent contamination and evaporation of the methylmethacrylate. Clean the thread of the bottle occasionally using a dry cloth to prevent the cap sticking.
- Do not use for long-term application, e.g. varnishing denture materials.
- Avoid liquid coming into contact with the skin. If there is accidental contact, wash with soap and water. If there is contact with the eyes, rinse immediately with plenty of water and consult a doctor.

Interactions

Avoid materials containing eugenol, moisture and oleaginous air. These can prevent polymerization at the contact point.

Contraindications

The application of *Varnish & Bond 2* is contraindicated when allergies against any components (especially methyl methacrylate) are known.

Side effects

Side effects are not known to date. In individual cases, hypersensitivity or contact allergy to components of the material cannot be excluded.

Technical data

Wavelength LED/halogen	approx. 450 nm
Curing time LED/Halogen ($\geq 400 \text{ mW/cm}^2$)	20 – 40 s
Wavelength UV – Light-curing unit	approx. 380 nm
Curing time UV – Light-curing unit	90 – 180 s

Composition

Multifunctional acrylates, methyl methacrylate, catalysts, stabilizers and additives.

Storage

Store at temperatures between $2^\circ\text{C} - 25^\circ\text{C} / 36^\circ\text{F} - 77^\circ\text{F}$. Please close bottle immediately after use to prevent volatilization of the methyl methacrylate. Do not use after date of expiry.

Available product sizes

Varnish & Bond 2

REF 2086 Bottle @ 5 ml

25 Brushes

Keep away from children! For dental use only!