»Tray adhesive for VPS«

de

Produktbeschreibung

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en-US

Product description

Tray adhesive for VPS is suitable for use on impression trays with minimal or no retention areas. It offers a secure adhesion of addition curing silicones to metal and plastic impression trays. When using a plastic impression tray, first check its compatibility with the tray adhesive.

Indication

All dental situations requiring final impressions e.g., tooth loss, caries and trauma.

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

Tray adhesive for VPS offers a secure adhesion of addition curing silicones to metal and plastic impression trays.

Limitations of use

Tray adhesive for VPS should only be used with addition curing, silicone-based impression materials.

Intended user

Dentist, dental assistant staff

Patient target group

People treated in the course of a dental procedure.

Recommended use

1. The impression tray should be clean, fat free and dry.
2. Use a brush to apply a thin, even layer of *Tray adhesive for VPS* to the inside of the impression tray.
3. Let the adhesive set for 3 min to ensure secure adhesion.
4. Add the impression material to the tray within 30 min to allow the adhesive to develop its full adhesive strength.

Notes for use

- Highly flammable
- Avoid breathing in the fumes
- Do not swallow

Interactions

When using a plastic impression tray, first check its compatibility with the tray adhesive.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!

- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Viscosity: 0.15 – 5 Pa s., secure bonding

Storage/disposal

- Store at room temperature (15 – 25 °C / 59 – 77 °F).
- Close bottles immediately after use to prevent contamination as well as evaporation of the solvent.
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Isoalkane, siloxane copolymer, modified resin, additives. Contains hexane

Packaging

REF 9046MDR	1 Bottle @ 10 ml
-------------	------------------