## »Bracket adhesive etch LC«

de

### Produktbeschreibung

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en-US

### Product description

Bracket adhesive etch LC is a bracket adhesive for use on etched enamel and consists of an adhesive paste, a primer and an etching gel. Curing is initiated by light.

### Indications

Malformed bite, malocclusion

### Contraindications

Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended use

Bracket adhesive for luting metal and ceramic brackets to etched enamel.

### Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.

### Intended users

Dentist

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Recommended use

#### Pretreat the enamel surface

1. Clean the enamel surfaces selected for luting the brackets with fluoride-free polishing paste and water.
2. Rinse teeth thoroughly with water; dry with oil- and moisture-free compressed air.

#### Apply the etching gel

3. Apply the etching gel to the enamel surface.
4. Let the etching gel take effect for 20 – 60 seconds.
5. Thoroughly rinse off the etching gel with water. Rinse each tooth for at least 5 seconds.
6. Dry tooth surfaces with oil-free and water-free air.
7. The etched enamel surface should appear chalky and opaque. If not, repeat the etching process. Protect the etched enamel surface from moisture and contamination.

#### Apply primer

8. Place 1 – 2 drops of primer in the hollow of a mixing pad.
9. Apply a thin layer of the primer to the etched enamel surface using a single-use brush.

Do not light-cure the primer. Light-curing takes place together with the bracket paste.

**Attention:** Too much primer can make attachment of the bracket difficult and have a negative effect on the adhesion values.

### Apply bracket adhesive paste

10. Slowly and evenly apply a small amount of bracket adhesive paste to the metal or ceramic bracket base.
11. Move the brackets to the desired final position, correct within 15 s if necessary and then press firmly in place.
12. Remove excess material carefully with a suitable instrument without changing the position of the bracket.

### Light-curing of the bracket adhesive

#### a) Metal bracket

- Light-cure the bracket adhesive twice for 10 s; once from above the metal bracket for 10 s and once from below the metal bracket for 10 s.
- Position the light emission tip of the curing light accordingly and hold it as close as possible (approx. 2 mm) to the upper or lower edge of the metal bracket.

#### b) Ceramic bracket

- Light-cure the bracket adhesive once for 10 s through the ceramic bracket. Hold the light tip approx. 2 mm away from the bracket and direct the light beam through the ceramic bracket.

### Remove residual bracket adhesive

13. After removing the bracket, carefully remove residual bracket adhesive using a scaler or a hard-metal finishing bur under dry working conditions, and using a continuous air stream.
14. Polish tooth surfaces using a rubber polisher and pumice powder.
15. Fluoridate tooth surfaces.

### Notes for use

- Avoid premature contamination of the adhesive paste and the primer.
- Light-curing units should emit at 450 nm and their light output should be checked regularly. The luminous intensity should be at least 1000 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the material.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

### Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration. Leftover methacrylate-based materials may influence the setting behavior of silicone materials.

### Residual risks/side effects

Adhesive paste and primer may trigger allergies. If there are known allergies to acrylates or methacrylates, do not use the material.

### Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.

## »Bracket adhesive etch LC«

- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Product performance characteristics

SBS bonding properties on etched enamel measured according to DIN 13990:2017 with metal or ceramic brackets are  $\geq 10$  MPa

### Storage/disposal

- Store in a dry place at room temperature (15 - 25 °C/59 - 77 °F)!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

### Composition

Etching gel: Phosphoric acid, water, SiO<sub>2</sub>, glycerol, additives.

Primer: Bis-GMA, TEDMA, polymer, CQ, additives.

Adhesive paste: Dental glass, bis-GMA, TEDMA, NaF, CQ, additives.

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

### Packaging

REF 1064MDR	1 Syringe Bracket adhesive paste, 5 ml 1 Bottle Primer, 5 ml 1 Bottle Etching gel, 10 ml 1 Spatula 50 Single-use Brushes 1 Mixing Pad
-------------	--