»Composite core build-up z dual«

de

Produktbeschreibung

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en-US

Product description

Composite core build-up z dual is a dual-cure, radiopaque core build-up composite that is automatically mixed and directly applied. Using the Intraoral Tips provided, even difficult regions are easily accessed, such as in root post cementation.

Indication

Severe loss of dental hard tissue

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

- All types of core build-ups and build-up fillings
- Cementation of root posts

Limitations of use

Do not use directly on the exposed pulp

Intended user

Dentist

Patient target group

People treated in the course of a dental procedure.

Recommended use

Cavity preparation

Note: It is advisable to use a rubber dam.

1. Before applying *Composite core build-up z dual*, remove caries and existing fillings and provide for appropriate pulp/dentin protection.
2. Irrigate the cavity carefully and air dry. Avoid both excessive residual moisture and overdrying. For adhesive anchorage of *Composite core build-up z dual* to the hard tissue, a dual-cure bonding agent should be used (e.g., *Total-etch adhesive dual*). Observe the instructions for use provided by the manufacturer.

Core build-up

1. To facilitate application of *Composite core build-up z dual*, place a matrix if necessary around the core that is to be prepared.
2. Apply *Composite core build-up z dual* directly inside the preparation. The material can be formed using commercially available instruments, e.g. a Heidemann spatula. The working time of *Composite core build-up z dual* is a maximum of 90 s. Curing can be accelerated by exposure to a suitable lamp. The light-curing time is 20 s at a light intensity of $\geq 1000 \text{ mW/cm}^2$ or 40 s at a light intensity of $\geq 600 \text{ mW/cm}^2$. With a chemical self-cure only, the restoration must be allowed to cure for at least 360 s before it can support a full load.

3. Remove the matrix after *Composite core build-up z dual* has cured. If a pre-fabricated core build-up form is used, *Composite core build-up z dual* can be applied directly into the core build-up form.
4. Then place the core build-up form onto the prepared stump and keep in position until *Composite core build-up z dual* has cured.

Cementation of root posts

1. Prepare the root post in accordance with the instructions for use provided by the manufacturer.
2. Following preparation of the root canal, etch the inside of the canal as per the instructions of the manufacturer of the bonding agent used.
3. Work in the bonding agent. Use a paper tip to remove any excess prior to light-curing.
4. Apply *Composite core build-up z dual* inside the prepared canal. The root post can be coated with *Composite core build-up z dual* before it is placed into the filled canal.
5. Once the root post has been positioned, light-cure the material with a suitable lamp for 20 s at a light intensity of $\geq 1000 \text{ mW/cm}^2$ or 40 s at a light intensity of $\geq 600 \text{ mW/cm}^2$ to fix the root post. Thanks to the dual curing of *Composite core build-up z dual*, the areas inaccessible to light cure chemically within 360 s. *Composite core build-up z dual* can also be used for the core build-up of the coronal portion of the tooth (see the section on "Core build-up"). A matrix should be used here due to the low viscosity of the material.

Fabricating a temporary restoration

If a resin-based temporary restoration is to be placed onto the prepared stump, use a suitable separation medium (e.g., Vaseline) for isolation.

Timing

Maximum working time including mixing time	90 s
Light-curing time with a light intensity of $\geq 1000 \text{ mW/cm}^2$	20 s
Light-curing time with a light intensity of $\geq 600 \text{ mW/cm}^2$	40 s
Minimum intraoral curing time including mixing time	360 s

Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light units should emit at 450 nm and their light output should be checked regularly. The luminous intensity should be at least 600 mW/cm^2 . Place the light-curing unit as close as possible to the material.
- The material should be at room temperature before use. Cooler material cures more slowly and has a higher viscosity.

Interactions

- The use of *Composite core build-up z dual* with single-bottle bonding systems is contraindicated.
- Do not use *Composite core build-up z dual* together with materials that contain eugenol!

»Composite core build-up z dual«

Using the double-barreled syringe

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new double-barreled syringe, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the double-barreled syringe as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment.

Note: When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

The product complies with ISO 4049:2019, Type 2, Class 3.

Storage/disposal

- Store at 2 – 25 °C/35.6 – 77 °F.
- Leave the used mixing tip on the cartridge or syringe as a cap after application!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Dental glass, silica (sil.), Bis-GMA, UDMA, HDDMA, DDDMA, MMHE, additives.
Inorganic fillers: approx. 40 vol% (0.02 – 2.5 µm).

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

Packaging

REF 2100MDR	1 Double-barreled syringe (A3) @ 5 ml
	10 Mixing tips
	10 Intraoral Tips