## »Flowable composite«

de

### Produktbeschreibung

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en-US

### Product description

*Flowable composite* is a light-curing, radiopaque restoration material for restorative treatment.

### Indications

Loss of tooth substance due to caries, trauma, tooth attrition/abrasion or development disorders.

### Contraindications

- Do not use directly on the exposed pulp.
- Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended use

As the first layer for class I and II fillings, anterior teeth fillings (class III, IV), class V fillings (cervical caries), small fillings of all types, for blocking out undercuts, for repairing composite/ceramic veneerings, for repairing bis-acrylate temporary restorations.

### Limitations of use

Application of *Flowable composite* is contraindicated if a dry working area or the recommended application technique is not possible.

### Intended users

Dentist

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Recommended use

#### Recommended use for fillings

1. Select the color of the composite to be used.
2. Provide for a dry operating field, preferably with a rubber dam.
3. Prepare the cavity in accordance with the rules for the adhesive technique.
4. Clean the cavity with water spray.
5. Dry the cavity with air.
6. Selectively apply pulp protection (calcium hydroxide) only if required (near the pulp) and cover with a pressure-resistant cement.
7. Place the matrix/interdental wedge.
8. Carry out conditioning and application of the bonding agent according to the instructions for use for the product used.
9. Apply the *Flowable composite* in layers of maximum 2 mm thickness.
10. Polymerize each layer for 40 seconds with a light intensity of 600 mW/cm<sup>2</sup> or for 20 seconds with a light intensity of 1000 mW/cm<sup>2</sup>. If a metal matrix is used, additionally expose the buccal or lingual/palatal surface after removal of the matrix.
11. Finish with suitable finishing-burs or fine-grained diamonds.
12. Perform an occlusion check.
13. Polish with silicone polishers as well as polishing discs and polishing strips.

### Recommendations for repairing temporary restorations

For the repair of an in situ temporary restoration, lightly roughen the site to be repaired using a bur or a sandblaster and apply a bonding agent. When repairing a freshly-made temporary restoration, it is not necessary to roughen the surface or use a bonding agent.

### Notes for use

- If the material is applied into the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be at least 600mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the filling material.
- If another layer of filling material is to be added, do not remove the oxygen-inhibited layer, as this will be required for bonding with the next layer.
- Follow the instructions for use for other materials that you use in conjunction with *Flowable composite*.

### Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Product performance characteristics

The product complies with ISO 4049:2019, Type 2, Class 2, Group 1; Polymerization depth  $\geq 2$  mm by using a light-curing time of 40 seconds with a light intensity of  $\geq 600$  mW/cm<sup>2</sup> or by using a light-curing time of 20 seconds with a light intensity of  $\geq 1000$  mW/cm<sup>2</sup>

### Storage/disposal

- Store at room temperature (15-25 °C / 59-77 °F).
- Close syringe immediately after use.
- Do not use after the expiration date.

## »Flowable composite«

- Disposal must comply with national regulations.

### Composition

Dental glass, bis-GMA,UDMA, HDDMA, DDDMA, SiO<sub>2</sub>, additives.

Inorganic filling materials: approx. 41 vol% (0.02 – -3 µm).

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

### Packaging

REF 2035MDR	2 Syringes (A2) @ 1.5 g 10 Luer-Lock-Tip
-------------	---