»Nano hybrid composite«

de

Produktbeschreibung

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en-US

Product description

Nano hybrid composite is a light-cured, radiopaque filling material for restorations of all cavity classes.

Indication

Loss of dental enamel due to caries, trauma, tooth abrasion or development disorders

Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

Restorations of all cavity classes

Limitations of use

Application of *Nano hybrid composite* is contraindicated if a dry working area or the recommended application technique are not possible.

Intended user

Dentist

Patient target group

People treated in the course of a dental procedure.

Recommended use

1. Clean the tooth. Remove all cleaning residue with water spray. Determine the VITA shade. The use of a rubber dam is recommended.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
3. If necessary, etch using your preferred etching technique such as Total-Etch. Observe the manufacturer's instructions for the bonding agent used.
4. Apply the bonding agent in accordance with the manufacturer's instructions.
5. If necessary, apply a base of low-viscosity composite or compomer. We recommend using a matrix.
6. Apply the desired shade of *Nano hybrid composite* to the cavity with a maximum layer thickness of 2 mm. Polymerize each layer for 40 s with a light intensity of at least 600 mW/cm² or for 20 s with a light intensity of at least 1000 mW/cm².
7. Removal of any excess material, contouring, and finishing are done using diamond finishing tools and flexible polishing disks of decreasing grain size. *Nano hybrid composite* can be given a high shine using polishing paste together with polishing cups or discs.

Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm². Place the light-curing unit as close as possible to the filling material.
- Do not remove the oxygen inhibited layer, because this is required for the bond with the next layer.
- Follow the manufacturer's instructions for any other materials that are used with *Nano hybrid composite*.

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

The product complies with ISO 4049:2019, Type 1, Class 2, Group 1.

Polymerization depth: ≤ 2 mm by using a light-curing time of 40 s with a light intensity of ≥ 600 mW/cm² or by using a light-curing time of 20 s with a light intensity of ≥ 1000 mW/cm²

Storage/disposal

- Store at room temperature (15 – 25 °C/59 – 77 °F).
- Protect against exposure to direct sunlight!
- Close syringe/composite tip immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Dental Glass, EBPADMA, bis-GMA, NPG2PODA, UDMA, SiO₂, additives.

Inorganic fillers: approx. 65 vol% (0.02 – 0.7 µm).

»Nano hybrid composite«

Packaging

REF 2054MDR	1 Syringe (A2) @ 3 g
REF 2061MDR	20 Composite Tips (A2) @ 0.3 g