## »Nano hybrid flowable composite«

de

### Produktbeschreibung

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en-US

### Product description

*Nano hybrid flowable composite* is a light-curing, radiopaque restoration material for restorative treatment.

### Indication

Loss of dental enamel due to caries, trauma, tooth abrasion or development disorders

### Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended purpose

As the first layer for class I and II fillings, anterior teeth fillings (class III, IV), class V fillings (cervical caries), small fillings of all types, blocking out undercuts, repairing composite/ceramic veneers, repairing bis-acrylate temporary restorations

### Limitations of use

Application of *Nano hybrid flowable composite* is contraindicated if a dry working area or the recommended application technique is not possible.

### Intended user

Dentist

### Patient target group

People treated in the course of a dental procedure.

### Recommended use

1. Determine the composite shade to be used.
2. Isolate the lesion to maintain a dry operating field. The use of a rubber dam is recommended.
3. Prepare the cavity in accordance with the rules for the adhesive technique.
4. Clean the cavity with water spray.
5. Dry the cavity with oil-free air.
6. Selectively apply pulp protection (calcium hydroxide) only if required (near the pulp) and cover with a pressure-resistant cement.
7. Place the matrix/interdental wedge.
8. Apply the bonding agent according to the manufacturer's instructions for the product being used.
9. Apply *Nano hybrid flowable composite* in layers of 2 mm thickness at maximum.
10. Polymerize each layer for 40 s with a light intensity of at least 600 mW/cm<sup>2</sup> or 20 s with a light intensity of at least 1000 mW/cm<sup>2</sup>. When using a metal matrix, additionally light cure the buccal or lingual/palatal regions after removing the matrix.
11. Complete the restoration with suitable finishing burs or fine-grained diamond polishers.
12. Perform an occlusion check.
13. Polish with silicone polishers or polishing discs and strips.

### Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm<sup>2</sup>. Place the light-curing unit as close as possible to the filling material.
- Do not remove the oxygen-inhibited layer if a further layer of filling material is to be applied, because this is required for the bond with the next layer.
- Follow the manufacturer's instructions for any other materials that are used with *Nano hybrid flowable composite*.

### Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Product performance characteristics

The product complies with ISO 4049:2019, Type 2, Class 2, Group 1. Polymerization depth  $\geq 2$  mm by using a light-curing time of 40 s with a light intensity of  $\geq 600$  mW/cm<sup>2</sup> or by using a light-curing time of 20 s with a light intensity of  $\geq 1000$  mW/cm<sup>2</sup>

### Storage/disposal

- Store in a dry place at room temperature (15 – 25 °C/59 – 77 °F).
- Close syringe immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.

### Composition

Dental glass, YbF<sub>3</sub>, EBPADMA, bis-GMA, TEDMA, UDMA, SiO<sub>2</sub>, additives.  
Inorganic fillers: approx. 38 vol% (0.02 – 0.7 µm).

## »Nano hybrid flowable composite«

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

### Packaging

REF 2060MDR	2 Syringes (A2) @ 1.5 g, 10 Luer-Lock Tips
-------------	---