»Self adhesive composite luting cement II«

de

Produktbeschreibung

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en-US

Product description

Self adhesive composite luting cement II is a universal, radiopaque, self-adhesive, dual-curing composite luting cement for cementing indirect restorations and root posts. When using *Self adhesive composite luting cement II*, there is no need to apply any additional bonding agent to the tooth substance.

Indication

- Extensive tooth defects: Tooth decay, severe wear, traumatic injuries of the dental enamel
- Tooth loss

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

- Permanent cementation of crowns, bridges, inlays and onlays made of metal, ceramics or composites.
- Permanent cementation of root posts made from fiber-reinforced composites.
- Permanent cementation of crowns and bridges on implant abutments comprised of titanium or zirconium oxide.
- Permanent cementation of electroformed caps and restorations made of titanium, PEEK or CoCrMo to tertiary structures or attachments, or to telescopic elements comprised of metals or methyl methacrylates.

Limitations of use

- *Self adhesive composite luting cement II* is contraindicated in cases where the preparation does not offer sufficient retention (e.g., veneers).
- Do not apply *Self adhesive composite luting cement II* onto the exposed pulp or dentin close to the pulp

Intended user

Dentist

Patient target group

People treated in the course of a dental procedure.

Recommended use

Application

1. Remove the temporary restoration and carefully clean off any residual temporary luting cement from the cavity or prepared core. Rinse with a water spray. Dry using oil-free air, but avoid overdrying the dentin.
2. Verify the accuracy of the fit of the restoration.
3. Then clean the restoration thoroughly and pretreat as follows:
 - **Etchable glass ceramics:** Etch the restoration with 5 % hydrofluoric acid in accordance with the manufacturer's instructions, rinse thoroughly using a water spray and dry using oil-free air, and apply silane (e.g., *Silane Primer*) in accordance with the manufacturer's instructions.

- **Oxide ceramics, metals and metal alloys:** Pretreat in accordance with the manufacturer's instructions. Unless otherwise stated, sandblast the inside of the restoration with aluminum oxide (< 50 µm), clean with alcohol, and dry with air that is free of oil and water.
 - **Fiber-reinforced root posts:** Clean the root post with ethanol, dry with oil-free air, and apply silane (e.g., *Silane Primer*) in accordance with the manufacturer's instructions.
 - **PEEK:** Pretreat in accordance with the manufacturer's instructions. If not described otherwise, clean the inner side of the restoration using isopropyl alcohol, sandblast (aluminum oxide, 50 – 110 µm, 2 – 3 bar), and clean with air free of oil and water. Apply *Varnish & Bond II* as specified by the manufacturer and light cure.
 - **Restorations made from other materials:** Pretreat in accordance with the manufacturer's instructions.
4. If required, clean the cavity or preparation thoroughly using a water spray to remove any residual temporary cement, and then dry. In doing so, ensure that the dentin retains a little residual moisture and does not become too dry.
 5. Apply *Self adhesive composite luting cement II* directly into the restoration using the tip. The working time of *Self adhesive composite luting cement II* is approx. 1 min. Place the restoration fully into position using slight pressure. Allow the patient to slowly adopt the habitual occlusion. The restoration can support a full load after 7 min.
 6. If necessary, light-cure the material with a suitable light-curing unit for at least 20 s at 1000 mW/cm² or 40 s at 600 mW/cm².

Removing excess material

Self-curing

The excess should be removed in the elastic phase, i.e. approx. 30 – 60 s after application, using a microbrush/brush/foam pellets/scaler or dental floss. It is particularly important to ensure that excess is removed quickly in areas that are difficult to access (proximal, gingival margins).

Self-curing with additional light-curing

Light cure the excess material briefly (1 – 2 s) so that it can be easily removed with a scaler. It is particularly important to ensure that excess is removed quickly in areas that are difficult to access (proximal, gingival margins), as it is difficult to remove excess from these areas once the material has fully cured.

Notes for use

- After it has been expressed from the double-barreled syringe, curing of *Self adhesive composite luting cement II* can be accelerated using strong lighting.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be at least 600 mW/cm². Place the curing light as close as possible to the material.

Interactions

Do not use *Self adhesive composite luting cement II* together with materials that contain eugenol. The use of temporary non-eugenol cements (e.g., non-eugenol *Temporary luting cement*) is recommended.

»Self adhesive composite luting cement II«

Using the double-barreled syringe

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new double-barreled syringe, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the double-barreled syringe as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment.

Note: When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

The product complies with DIN ISO/TS 16506:2021 Class 3 (including shade stability).

Storage/disposal

- Store at 2 – 8 °C/36 – 46 °F in a cool, dry place!
- Protect against exposure to direct sunlight!
- After application, leave the used mixing tip on the double-barreled syringe as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

Dental glass, EBPADMA, UDMA, TEDMA, HEMA, MDP, SiO₂, bis-GMA, additives.

Filling materials: approx. 50 vol% (0.02 – 3.0 µm).

Packaging

REF 2098MDR	1 Double-barreled syringe @ 5 ml, shade Universal 10 Mixing tips
-------------	--