»Bis-acrylic temporary crown & bridge material 1:1«

de

Produktbeschreibung

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en-us

Product description

Bis-acrylic temporary crown & bridge material 1:1 is a methacrylate-based temporary crown and bridge material. The material is automatically mixed in a ratio of 1:1. Thanks to its excellent mechanical values, *Bis-acrylic temporary crown & bridge material 1:1* is extremely fracture-resistant and thus most suitable for fabricating long-term temporary restorations. *Bis-acrylic temporary crown & bridge material 1:1* guarantees a highly accurate fit of the finished temporary restoration with optimal marginal fit.

Indication

- Extensive tooth defects: Tooth decay, severe wear, traumatic injuries of the dental enamel
- Shape and positional anomalies of the teeth
- Tooth loss.

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns, veneers and long-term provisional restorations

Limitations of use

Use the material to create temporary restorations only

Intended user

Dentist, dental assistant staff and dental technician

Patient target group

People treated in the course of a dental procedure.

Recommended use

A dental impression should be taken using alginate or silicone before preparing a stump for a crown or bridge or before a planned extraction. The interdental flash of the impression should be cut away to improve the stability of the temporary restoration to be fabricated at a later stage. When using alginate impressions, apply the *Bis-acrylic temporary crown & bridge material 1:1* as soon as possible after taking the impression.

Recommended use and elastic phase

1. Dry the prepared teeth and lightly coat the stumps, surrounding tissue and any existing synthetic material with a separating agent (e.g., Vaseline).
2. Apply *Bis-acrylic temporary crown & bridge material 1:1* into the deepest section of the impression and then fill up to the gingival areas. Always keep the mixing tip immersed in the material to avoid air bubbles.
3. Reposition the filled impression in the mouth within 50 s.
4. *Bis-acrylic temporary crown & bridge material 1:1* has a firm, elastic consistency approx. 105 – 135 s after mixing and can be easily removed from the patient's mouth. The curing process must be monitored intra-orally (e.g., with a probe), as the temporary restoration can only be easily removed during the elastic phase.

5. The final chemical curing of *Bis-acrylic temporary crown & bridge material 1:1* occurs within 300 s after mixing. To this end, the temporary restoration should be placed back into the impression. After the removal of the oxygen-inhibited layer with a solvent (e.g., alcohol), the temporary restoration can be refined and polished. For a shiny surface sealing, *Varnish & Bond 2* can be applied to the entire temporary restoration and light-cured (according to the manufacturer's instructions).
6. Commercially available temporary cements (e.g., *Temporary luting cement*) can be used for luting the temporary restoration.

Repair of temporary restorations

A feature of temporary restorations fabricated using *Bis-acrylic temporary crown & bridge material 1:1* is their high mechanical strength. If, however, a temporary restoration should fracture, we recommend the following procedure:

Fracture of the temporary restoration shortly after fabrication

Bond the fractures with freshly mixed *Bis-acrylic temporary crown & bridge material 1:1*.

Press the fragments together for approx. 180 s. After *Bis-acrylic temporary crown & bridge material 1:1* is completely chemically cured, the temporary restoration can be processed as usual.

Fracture of an older or in situ temporary restoration

Roughen the fractured sections and prepare undercuts.

Bond the prepared fracture with freshly mixed *Bis-acrylic temporary crown & bridge material 1:1*. The use of a suitable bonding agent (e.g., *Varnish & Bond 2*) for improved adhesion is recommended.

Press the fragments together for approx. 180 s.

After *Bis-acrylic temporary crown & bridge material 1:1* is completely chemically cured, the temporary restoration can be processed as usual.

Timing

Insert the impression into the mouth	0 – 50 s
Removal from the mouth	105 – 135 s
Removal of the oxygen-inhibited layer with solvent	≥ 300 s

Notes for use

If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Note: The material should be at room temperature before use. Cooler material cures more slowly and has a higher viscosity.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of silicone impression and bite registration materials.

»Bis-acrylic temporary crown & bridge material 1:1«

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Cartridge handling

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new cartridge, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the cartridge as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment.

Note: When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Compressive strength: ≥ 350 MPa, flexural strength ≥ 100 MPa

Storage/disposal

- Store at 2 – 25 °C /36 – 77 °F in a dry place!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Dental glasses, UDMA, EBPADMA, bis-GMA, TEDMA, additives.

Packaging

REF 1108MDR	1 Cartridge (A2) @ 50 ml 10 Mixing tips
-------------	--