»Bis-acrylic temporary crown & bridge material 4:1«

de

Produktbeschreibung

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Dentalglas, EBPADMA, Urethanacrylat oligomer, Polyesterharz, Urethanmethacrylat, SiO₂, Additive.

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en-US

Product description

Bis-acrylic temporary crown & bridge material 4:1 is a bis-acrylic based temporary crown and bridge material that is free from methyl methacrylate. *Bis-acrylic temporary crown & bridge material 4:1* offers special "Snap-Set" setting behavior and consequently has a long plastic working stage, during which the temporary restoration can be easily removed from the patient's mouth. *Bis-acrylic temporary crown & bridge material 4:1* is exceptionally fracture-resistant and color stable due to its outstanding mechanical properties. and is therefore ideal for the fabrication of long-term provisional restorations. *Bis-acrylic temporary crown & bridge material 4:1* guarantees a highly accurate fit of the finished temporary restoration with optimal marginal integrity. *Bis-acrylic temporary crown & bridge material 4:1* is mixed automatically.

Indication

- Extensive tooth defects: Tooth decay, severe wear, traumatic injuries of the dental enamel
- Shape and positional anomalies of the teeth
- Tooth loss

Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

Intended purpose

Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns, veneers and long-term provisional restorations.

Limitations of use

Use the material to create temporary restorations only.

Intended user

Dentist, dental assistant staff and dental technician

Patient target group

People treated in the course of a dental procedure.

Recommended use

A dental impression should be taken using alginate, silicone or a thermoplastic impression material before preparing for a crown or bridge or before a planned extraction. A prepared thermoformed splint or a strip crown may also be used as a mold. The interdental flash of the impression should be cut away to improve the stability of the temporary restoration to be fabricated at a later stage. When using alginate impressions, use the *Bis-acrylic temporary crown & bridge material 4:1* as soon as possible after taking the impression.

Recommended method for extraoral light-curing

1. Dry the prepared teeth and lightly coat the stumps, surrounding tissue and any existing synthetic material with a separating agent (e.g., Vaseline).

2. Apply *Bis-acrylic temporary crown & bridge material 4:1* into the deepest section of the impression and then fill up to the gingival areas. Always keep the mixing tip immersed in the material to avoid air bubbles.
3. Reposition the filled impression in the mouth within 40 s.
4. *Bis-acrylic temporary crown & bridge material 1:1* has a firm, elastic consistency approximately 90 – 140 s after mixing and can be easily removed from the patient's mouth. The curing process must be monitored intra-orally (e.g., with a probe), as the temporary restoration can only be easily removed during the elastic phase.
5. The material is fully chemically cured after 300 s. *Bis-acrylic temporary crown & bridge material 4:1* can be light-cured to accelerate curing. For light-curing, the material should be exposed for at least 20 s per crown or pontic at an exposure intensity of ≥ 1000 mW/cm² or 40 s at an exposure intensity of ≥ 600 mW/cm² using a hand lamp. For more extensive work, the entire temporary restoration should be exposed in a dental polymerization unit for at least 60 s.
6. The oxygen inhibition layer can be removed with a solvent (e.g., alcohol) after curing. The temporary restoration can be finished and polished. For a high quality final result, the temporary restoration can also be treated with a high-gloss surface sealer (e.g., *Varnish & Bond*).
7. Commercially available temporary cements (e.g., *Temporary luting cement*) can be used for luting the temporary restoration.

Recommended use for intraoral light-curing (veneer)

Important: Intraoral light-curing should only be performed with veneers made with *Bis-acrylic temporary crown & bridge material 4:1*. These should have a maximum thickness of 1 mm. All other temporary restorations must be cured extraorally!

A transparent impression material (silicon impression material, thermoplastic impression polymer, thermoforming film or strip crown) is required for intraoral light-curing.

1. Dry the prepared teeth and lightly coat the stumps, surrounding tissue and any existing synthetic material with a separating agent (e.g., Vaseline).
2. Apply *Bis-acrylic temporary crown & bridge material 4:1* into the deepest section of the impression and then fill up to the gingival areas. Always keep the mixing tip immersed in the material to avoid air bubbles.
3. Reposition the filled impression in the mouth within 40 s.
4. Expose the material directly through the impression for at least 20 s at an exposure intensity of ≥ 1000 mW/cm² or 40 s at an exposure intensity of ≥ 600 mW/cm² per unit using a hand lamp. When doing so, please keep in mind that the distance of the light-curing unit has an important influence on curing speed. Light-curing time must be increased if necessary.
5. When *Bis-acrylic temporary crown & bridge material 4:1* is fully cured, remove the temporary restoration from the mouth, remove the oxygen inhibition layer using alcohol and then finish and polish the temporary restoration extraorally.

»Bis-acrylic temporary crown & bridge material 4:1«

For a high quality final result, the temporary restoration can also be treated with a high-gloss surface sealer (e.g., *Varnish & Bond*).

- Commercially available temporary cements (e.g., *Temporary luting cement*) can be used for luting the temporary restoration.

Recommended use for intraoral light-curing with spot-etching technique

Important: Intraoral light-curing should only be performed with veneers made with *Bis-acrylic temporary crown & bridge material 4:1*. These should have a maximum thickness of 1 mm. All other temporary restorations must be cured extraorally!

A transparent impression material (silicon impression material, thermoplastic impression polymer, thermoforming film or strip crown) is required for intraoral light-curing.

- Dry the prepared teeth and use the spot-etching technique.
- Apply *Bis-acrylic temporary crown & bridge material 4:1* into the deepest section of the impression and then fill up to the gingival areas. Always keep the mixing tip immersed in the material to avoid air bubbles.
- Reposition the filled impression in the mouth within 40 s.
- Expose the material directly through the impression for at least 20 s at an exposure intensity of $\geq 1000 \text{ mW/cm}^2$ or 40 s at an exposure intensity of $\geq 600 \text{ mW/cm}^2$ per unit using a hand lamp. When doing so, please keep in mind that the distance of the light-curing unit has an important influence on curing speed. Light-curing time must be increased if necessary.
- The impression can be removed from the mouth once the *Bis-acrylic temporary crown & bridge material 4:1* is fully cured. The oxygen inhibition layer can be removed using alcohol and the temporary restoration can then be finished and polished intraorally. For a high quality final result, the temporary restoration can also be treated with a high-gloss surface sealer suitable for intraoral use.

Repair of temporary restorations

A feature of temporary restorations fabricated using *Bis-acrylic temporary crown & bridge material 4:1* is their high mechanical strength. If, however, a temporary restoration should fracture, we recommend the following procedure: Fracture of the temporary restoration shortly after fabrication:

- Bond the fractures with freshly mixed *Bis-acrylic temporary crown & bridge material 4:1*.
- Press the fragments together and cure the freshly added material for at least 20 s at an exposure intensity of $\geq 1000 \text{ mW/cm}^2$ or 40 s at an exposure intensity of $\geq 600 \text{ mW/cm}^2$ using a hand lamp.
- Then finish and polish the temporary restoration as normal.

Fracture of an older or in situ temporary restoration:

- Roughen the fractured sections and prepare undercuts.
- Bond the prepared fracture with freshly mixed *Bis-acrylic temporary crown & bridge material 4:1*. We recommend using a composite bonding agent (e.g., *Varnish & Bond*) to improve the adhesion.

- Press the fragments together and cure the freshly added material for at least 20 s at an exposure intensity of $\geq 1000 \text{ mW/cm}^2$ or 40 s at an exposure intensity of $\geq 600 \text{ mW/cm}^2$ using a hand lamp.
- Then finish and polish the temporary restoration as normal.

Timing

Using only chemical self-curing	
Insert the impression into the mouth*	0 – 40 s
Removal from the mouth	90 – 140 s
Chemical self-curing is complete	$\geq 300 \text{ s}$

* Room temperature 23 °C including mixing time

Dual-curing	
Insert the impression into the mouth*	0 – 40 s
Removal from the mouth	90 – 140 s
Light-curing time with a light intensity of $\geq 600 \text{ mW/cm}^2$	40 s
Light-curing time per unit with a light intensity of $\geq 1000 \text{ mW/cm}^2$	20 s
Light-curing time of the entire restoration using a light oven	$\geq 60 \text{ s}$

* Room temperature 23 °C including mixing time

Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light units should emit at 450 nm and their light output should be checked regularly. The light intensity should be a minimum of 600 mW/cm^2 and the light should be positioned as closely to the material as possible.
- The times given in the working sequence for *Bis-acrylic temporary crown & bridge material 4:1* are valid at a room temperature of 23 °C or an oral temperature of 35 – 37 °C and a relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Cartridge handling

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new cartridge, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the cartridge as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment.

Note: When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

»Bis-acrylic temporary crown & bridge material 4:1«

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Compressive strength: ≥ 250 MPa dual-curing, Flexural strength: ≥ 75 MPa dual-curing

Storage/disposal

- Store at room temperature (15 – 25 °C / 59 – 77 °F).
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Dental glass, EBPADMA, urethane acrylate oligomer, polyester resin, urethane methacrylate, SiO₂, additives.

Inorganic filler content: 46 wt.% = 26 vol.%

The inorganic filler particles vary in size from 0.02 to 1.5 µm.

Packaging

REF 1102MDR	1 Cartridge @ 50 ml, 15 Mixing tips
-------------	--