## »Temporary composite luting cement transparent«

de

### Produktbeschreibung

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en-US

### Product description

*Temporary composite luting cement transparent* is a transparent, dual-cure composite luting cement that is indicated for temporary luting of temporary crowns, bridges, veneers, inlays, and onlays. *Temporary composite luting cement transparent* is automatically mixed in the double-barreled syringe and applied directly.

### Indication

- Extensive tooth defects: Tooth decay, severe wear, traumatic injuries of the dental enamel
- Tooth loss

### Contraindication

Do not use the material where there are contact allergies or known allergies to any of the components.

### Intended purpose

- Temporary luting of temporary crowns, bridges, inlays and onlays and temporary veneers.
- Temporary luting of implant-supported temporary prosthesis.

### Limitations of use

Do not apply the material directly to the exposed dental pulp.

### Intended user

Dentist

### Patient target group

People treated in the course of a dental procedure.

### Recommended use

1. If required, apply a separating agent to the tooth core, core build-up, or temporary (see also Notes for use).
2. Apply the material to the temporary restoration.
3. Insert the restoration within 60 s after mixing has begun and attach to the dry, prepared tooth using slight pressure.
4. If necessary, instruct the patient to adopt the habitual occlusion.
5. Remove cement residues within 75 – 120 s after inserting the restoration with the aid of a microbrush, brush, foam pellet, or scaler. Drying with oil-free and water-free air makes excess cement easier to identify. Remove interdental excess carefully using dental floss. As an alternative, it is possible to remove excess cement by briefly curing it using light (1 – 2 s per side) and then detaching it with a scaler.
6. Allow the material to chemically cure completely for 360 – 420 s after mixing has begun.
7. Translucent restorations: If necessary, additionally cure the material with a suitable light-curing unit for at least 20 s at a light intensity of at least 600 mW/cm<sup>2</sup> per side.

### Timing

Working time (including mixing time)	max. 60 s
Removal of excess material (without light-curing).	75 – 120 s after the restoration is inserted
Optional: Light-curing time for removal of excess material (per side)	1 – 2 s
Intraoral curing	360 – 420 s
Light-curing (per side)	20 s with a light intensity of at least 600mW/ cm <sup>2</sup>

### Notes for use

- Light units should have a light intensity of 600 mW/cm<sup>2</sup>, emit at 450 nm, and their performance should be checked regularly. The light-curing unit should be positioned as closely to the material as possible. Observe the manufacturer's instructions.
- The material adheres to composites. Apply a thin layer of a separating agent (e.g., Vaseline, K-Y Jelly) to existing composite-based fillings, bases or core buildups before using the material. In the case of further use of the temporary, apply a separating agent in advance.
- Ultrasound is helpful for the removal of the cement from the restoration.

### Interactions

Not to be used in conjunction with materials that contain eugenol/clove oil as these inhibit the hardening of the materials.

### Using the double-barreled syringe

The double-barreled syringes and cartridges enable the automatic dispensing and mixing, which allows a time-saving direct application. When using a new double-barreled syringe, remove the transport cap first. The mixing tip can then be fitted. After use and until the next treatment, leave the mixing tip on the double-barreled syringe as a seal. The used mixing tip should be replaced with a new tip when the device is used for the next treatment.

**Note:** When using a double-barreled syringe or cartridge for the first time, discard a small amount of the initially ejected material (e.g., a pea-sized amount).

### Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

### Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

## »Temporary composite luting cement transparent«

- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

### Product performance characteristics

Compressive strength  $\leq$  150 MPa

### Storage/disposal

- Store in a dry place at room temperature (15 – 25 °C/59 – 77 °F)!
- Protect against exposure to direct sunlight!
- After application, leave the used mixing tip on the syringe as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

### Composition

Dental glass, EBPADMA, aliphatic urethane methacrylate, unsaturated polyester resin blend, TEDMA, HEMA, SiO<sub>2</sub>, additives.

### Packaging

REF 2092MDR	1 Double-barreled syringe @ 5 ml, 10 Mixing tips
-------------	---